



WELCOME

Welcome to Revive Light Therapy® LipCare, light therapy system for treatment lip creases & wrinkles, temporary increase in circulation, and lip plumping.

SYSTEM INCLUDES

- Revive Light Therapy® Lip Care Device
- Removable Mouth Tray
- USB Charging Cord
- User Guide
- Storage Bag

Revive Light Therapy® LIP CARE DEVICE - LIGHT THERAPY SYSTEM FOR THE TREATMENT OF LIP CREASES & WRINKLES, TEMPORARY INCREASE IN CIRCULATION & LIP PLUMPING

CONTRAINDICATIONS

- Do not use this device if your skin is light sensitive you are currently using medication that may cause skin sensitivity to light. Doing so may result in photosensitive reaction. Discontinue use immediately if you experience any sensitivity to using this device.
- Contact your doctor prior to using this device if you are or suspect you are taking drugs known to cause light sensitivity.

WARNINGS

- Do not operate device where it could come into contact with water. Doing so could cause an electric shock or result in injury.
- Do not use or charge device with USB power adapter that is damaged or faulty. Doing so could cause in an electric shock, result in injury or damage the device.
- Do not use or charge device with damaged USB charging cord. Doing so could cause in an electric shock, result in injury or damage the device. Please contact customer service.
- Use only the supplied universal USB charging cord. Use of incorrect charging cord could damage the device.
- Do not attempt to take the device apart. Device is not serviceable by user. This will void the warranty.
- Stop using the device immediately if you think it is defective. Please contact customer service.
- Discontinue use if device overheats or becomes uncomfortably hot. Please contact customer service if you experience overheating of the device or if the device becomes uncomfortably hot to your skin.
- Do not store the device in direct sunlight or location with high temperatures as this may cause damage to the plastic coating.

PRECAUTIONS

- Always unplug device from wall outlet after charging is complete.
- Do not attempt to open or repair the device or power supply.
- No modification of this equipment is allowed.
- Do not use this device on infants and young children without physician approval.
- Contact your doctor if you have any medical questions regarding use of this device.
- Ask you doctor prior to using if you are pregnant or suspect to be pregnant or have any suspicious of cancerous lesions.

CUSTOMER SERVICE

You can reach customer service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time. Service can also be reached by email at info@aduroled.com.



DISCLAIMER

This product is not intended to diagnose, cure or prevent any disease. LED technologies Inc makes no claims, representations or warranties regarding the ability of this product to cure any physical, skin or mental condition. A qualified health professional should be consulted with regards to any condition requiring medical attention.

DEVICE OVERVIEW

Revive Light Therapy® Lip Care is a light emitting diode (LED) device that harnesses energy to plump lips and treat wrinkles around the lip area. The device incorporates natural light energy in the amber (650nm), Red (630nm), Deep red (660nm) and infrared (880nm) spectrums that are proven to reduce wrinkles, result in increased circulation and lip plumping.

HOW IT WORKS

Youthful skin is rich in collagen which gives the skin lift and support. Collagen works deep within the skin like columns of tissue that support the skin's surface. Over time skin produces less collagen and with exposure to sunlight healthy collagen becomes damaged. This can result in wrinkles specially around the lip area. The lights from this device penetrate deep within the layers of your skin diminishing creases and wrinkles. Easy to use, Revive Light Therapy® Lip Care is programmed a 3 minute auto shutoff function for professional anti-aging results with short simple daily treatments in the comfort of your home.

PRODUCT REPRESENTATION

Led Technologies Inc is constantly improving our products with the latest technology. As a result the actual product may differ slightly from the product shown in images on the packaging and in this user guide.

ASSEMBLY INSTRUCTIONS

1. Unpack your Revive Light Therapy® Lip Care
2. Separate lip device from mouth tray
3. Plug USB charging cord into connector on the device
4. Plug USB charging cord into a USB power adapter (not supplied) and into outlet for charging
5. Charge your Revive Light Therapy® Lip Care for 2-4 hours. Device is fully charged when all three indicator lights are solid.
6. Place device in protective storage bag when not in use
7. A fully charged device yields 30 minutes of use or 10 treatments.

LIGHT SENSITIVITY TEST

Test your skin for light sensitivity prior to initial treatment

1. Power on device
2. Place device on or within 1/4" of the skin of your forearm for 3 minutes
3. After 3 minutes check the area. If the area has turned red /redness and last for more than 2 hours your skin is light sensitive
4. Discontinue use and contact Customer Service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time

OPERATING INSTRUCTIONS

1. Fully Charged device
2. Wash and dry area prior to use
3. Place silicone-tray into mouth securing device into place
4. Power up the device
5. Leave device over treatment area for 3 minutes. The device will auto shut off after completion of 3 minutes treatment.
6. Place device into storage bag

HELPFUL TIPS FOR BEST RESULT

- Use daily for best result
- Do not use device after application of creams that contain SPF, as this will block light penetration
- Wipe device clean after each use

MAINTENANCE

- This is a medical device. It is important to keep it clean.
- This device has no serviceable parts
- Device only includes an ON/OFF button. All settings are set at the factory
- Charge device only with supplied USB charging cord
- Power indicator lights will begin to flash when battery levels are low, indicating it is appropriate to begin charging

CLEANING

1. Turn the device off and unplug from any power source
2. Pull apart device
3. Wash mouthpiece in lukewarm water and dry
4. Wipe the device surface with alcohol wipes or damp cloth to remove any residue from skin oils
5. Clean device body by using damp vlot or microfiber cloth. Microfiber storage bag can be used for cleaning and polishing the body of the device
6. Allow device to dry completely prior to use

TROUBLESHOOTING

- System is not turning on

Is the device full charged? (Power indicator lights will be on and solid)

If device still does not turn on contact Customer Service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time

- System is not charging

Is USB charging cord firmly inserted into the device? (Do Power indicator lights turn on?)

Try different USB power adaptor to ensure it is working properly

Try different wall-outlet to determine if there is power to the outlet

If device still does not charge contact Customer Service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time

- System shuts off abruptly, prior to completion 3-minutes treatment

System is designed to shut off automatically after 3 minutes

System may shut off automatically after multiple repeated uses if the device exceeds desired temperature. Allow device to cool for several minutes and try again

If device continues to shut off prior to 3 minutes please contact Customer Service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time

- Some lights appear turned off

The infrared (IR) light spectrum is not visible to the human eye. These lights will look as if they are not on but are operating properly. If the red visible lights are on then the device is operating as designed.

LIMITED WARRANTY INFORMATION

Congratulations on your purchase. We value you as a client and strive for the highest quality products earn your satisfaction. Please use this as your reference for any customer service or warranty needs you may have.

WHAT IS COVERED BY THE LIMITED WARRANTY?

This 6 month limited warranty covers defects in materials and workmanship in your Revive Light Therapy® Lip Care if you purchased it from an authorised partner or directly from the distributor. Your 6 month limited warranty is valid for you alone, the original purchaser, and cannot be transferred. If the product was a gift, you can ask the gift giver to process the warranty claim on your behalf.

WHAT IS NOT COVERED BY THE LIMITED WARRANTY?

This limited warranty does not cover damage, defects or failure caused by abuse or misuse, improper care, impact or insertion of foreign objects, use with a voltage connector or loss of device.

HOW LONG DOES THE LIMITED WARRANTY LAST?

The limited warranty begins from the date of purchase and expires 6 months later. The warranty period is not extended or replaced by Revive during the warranty period.

WHAT WILL Revive DO IF YOUR DEVICE IS DEFECTIVE?

If you discover a covered defect in your Revive Light Therapy® Lip Care before the 6 month warranty period expires contact Customer Service at +31 20 6463693 Monday - Friday 09:00 AM- 05:00 PM European Time. All limited warranty services must be authorised and performed by Revive.

WARRANTY LIMITATIONS & DISCLAIMERS

This warranty gives you specific legal rights and is in addition to your statutory rights. Your remedies for the breach of this limited warranty are limited solely to the repair and replacement of your device as set forth in this warranty statement. To the maximum extent permitted by applicable law. Revive disclaims all implied warranties. If applicable law does not permit the disclaimer of implied warranties then the duration of implied warranties are hereby expressly limited to the duration of the express limited warranty set forth above.

To the maximum extent permitted by applicable law in no event shall Revive be liable for any special, incidental, consequential, punitive, or indirect damages based upon breach of warranty, breach of contract, negligence, strict product liability former other legal theory even if Revive is advised or aware that such damages might arise. Such damages include, but are not limited to, loss of device, cost of any substitute equipment, the claims of third parties, or damage to property. If the applicable law does not allow limits on warranties or on remedies for breach in certain transactions, the limits in this paragraph and the disclaimer or exceptions from warranty may not apply. This limitation does not apply to claims for death or personal injury caused by a defect in the device or any other liability that cannot be excluded or limited by law.

These warranty terms, conditions, and limitations allocate the risk of product failure between Revive and you, the purchaser of Revive Light Therapy® Lip Care. This allocation is recognised by both parties and is reflected in the price of Revive Light Therapy® Lip care. By using Revive Light Therapy® Lip care you acknowledge that you have read these terms, conditions, and limitations, understand them and are bound by them.

Manufactured by



Greenwood Village, CO 80111
800.337.9565 Made in China

Our products are continually
improving. The actual product may
differ from the product shown.